TRANSPORT MEDIUM

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TRANSPORT MEDIUM

TM001/TM002/TM003/TM004/TM005/TM006/TM008/TM010/TM011/TM012/TM014/TM016/TM017/TM018/TM019/TM020/TM021/TM022: The transport medium avoids the drying up of the sample, improves microorganisms viability between collection and inoculation, and decreases bacterial overgrowth in contaminated samples.

WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use only.

The transport medium should appear as a transparent orange-coloured solution. If turbidity or a shift in colour to yellow or dark pink is observed, the medium must be discarded.

Check that the wrapping of the swab is not deteriorated to guarantee its sterility.

Once the sample has been collected, handle the material as potentially infectious.

Any serious incident that occurs in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

MATERIALS PROVIDED

VIRCELL TRANSPORT MEDIUM: 50 x 1 ml (REF: TM017) or 50 x 2 ml (REF: TM001/TM003/TM005/TM006/TM010/TM011/TM012/TM014/TM018/TM019/TM020) or 50 x 3 ml (REF: TM002/TM008/TM016/TM021/TM022) or 1 x 100 ml (REF: TM004) of Hank's BSS containing HEPES, gelatine, BSA, saccharose and compatible antibiotics.

COMMERCIAL SWABS: 100 x sterile Rayon swab (REF: TM001), 50 x sterile Rayon swab (REF: TM011), 50 x sterile Dacron swab (REF: TM003), 50 x sterile nylon microfiber swab (REF: TM008/TM012/TM014/TM016/TM017), 100 x sterile Dacron swab (REF: TM006), 50 x sterile Rayon swab + 50 x sterile nylon microfiber swab (REF: TM010/TM022), 50 x sterile nylon microfiber pediatric swab (REF: TM018/TM019) or 50 x sterile nylon microfiber swab + 50 x sterile nylon microfiber pediatric swab (REF: TM020).

STORAGE AND HANDLING CONDITIONS

Store at 2-25°C. Do not use the kit reagents beyond the expiration date. This will be valid only if reagents are stored closed and at 2-25°C.

IN-USE STABILITY

VIRCELL TRANSPORT MEDIUM: Refer to package label for expiration date (at 2-25°C).

VIRCELL, S.L. does not accept responsibility for the mishandling of the reagents included in the kit.

CONDITIONS FOR COLLECTION, HANDLING AND PREPARATION OF THE SPECIMEN

Samples for isolation must be collected by skilled personnel.

Collect samples during the acute phase of the disease.

Collect samples following established standard procedures.

It is very important to collect epithelial cells from the affected areas, hence samples must be collected by thoroughly swabbing or scraping these areas.

ASSAY PROCEDURE

Collect the sample with the sterile swab.

Put the swab into the tube containing the transport medium.

Cut off the end of the swab stick standing out of the tube, taking care not to touch the rim of the tube.

Close the tube tightly.

Fill in the data of the patient, type of sample and date and hour of collection.

Send to the laboratory for analysis at room temperature. For prolonged transport, use a refrigerated delivery service at 2-8°C.

If samples cannot be processed immediately upon receipt in the laboratory, they should be stored at 2-8°C up to 48 hours. For longer storage before inoculation, keep below -70°C. Avoid storage at -20°C since it has a deleterious effect on viability.

INTERNAL QUALITY CONTROL

Each batch is subjected to internal quality control (Q.C.) testing before batch release. Final Q.C. results for each particular lot are available.

LIMITATIONS OF USE

This product has been validated only in microorganisms culture techniques in laboratory. The use of the product in different applications must be previously validated by the end user.

PERFORMANCE CHARACTERISTICS

Viability assays were performed with *Chlamydomyphila pneumoniae*, *Chlamydia trachomatis*, enterovirus, cytomegalovirus and *Mycoplasma pneumoniae* in VIRCELL TRANSPORT MEDIUM against another available commercial the transport medium. 5 replicas of 3 different concentrations of each microorganism were assayed at 0, 24, 48, 72 h and 1 week, to 2-8°C and at room temperature. Viability was determined by counting fluorescent focus (CMV, *C. pneumoniae*, *C. trachomatis*, enterovirus) or colonies (*M. pneumoniae*). The assay result is showed in charts in the annex "PERFORMANCE".

Note: In the charts are represented only one of the dilutions used. Results of the other dilutions can be extrapolated.

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COMMERCIAL SWABS: 100 x tampone in Rayon sterile (REF: TM001), 50 x tampone in Rayon sterile (REF: TM011), 50 x tampone in Dacron sterile (REF: TM003), 50 x tampone in microfibra di nylon sterile (REF: TM008/TM012/TM014/TM016/TM017), 100 x tampone in Dacron sterile (REF: TM006), 50 x tampone in Rayon sterile + 50 x tampone in microfibra di nylon sterile (REF: TM010/TM022), 50 x tampone pediatrico in microfibra di nylon sterile (REF:TM018/TM019) o 50 x tampone in microfibra di nylon sterile + 50 x tampone pediatrico in microfibra di nylon sterile (REF: TM020).

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