

# ELISA KITS

## For antibody detection of infectious diseases

With more than 20 years of experience, Vircell offers ELISA kits for detecting specific antibodies against Herpes virus, TORCH and a wide range of infectious agents causing respiratory diseases and sexually transmitted diseases, among others, in human serum/plasma samples.

### MAIN BENEFITS

- Ready-to-use reagents.
- No sample predilution required.
- Colour-coded plates with individual break-apart wells.
- Semiquantitative protocols available for some parameters.
- $\mu$ -Capture technique for TORCH, Dengue and Parvovirus IgM determinations.
- G/M kits allow for customized testing of IgG and/or IgM.

### A SOLUTION FOR ALL

- Bacterial Respiratory Infections
- Viral Respiratory Infections
- TORCH Infections
- Vaccine-preventable Infections
- Vector-borne Infections
- Sexually Transmitted Infections
- Zoonotic Infections

### OTHER THINGS THAT MATTER

- Suitable for all automated ELISA systems.
- In-house production of high-quality antigens.
- Featured in more than 120 scientific publications.
- Vircell takes part in External Quality Control Programs.



# ELISA KITS FOR A WIDE RANGE OF INFECTIOUS DISEASES

## BACTERIAL RESPIRATORY INFECTIONS

<b>G/M1003</b>	CHLAMYDIA ELISA IgG/IgM
<b>A1007</b>	CHLAMYDOPHILA PNEUMONIAE ELISA IgA
<b>G1007</b>	CHLAMYDOPHILA PNEUMONIAE ELISA IgG (incl. calibrators)
<b>M1007</b>	CHLAMYDOPHILA PNEUMONIAE ELISA IgM
<b>G1001</b>	COXIELLA BURNETII ELISA IgG
<b>M1001</b>	COXIELLA BURNETII ELISA IgM
<b>G1000</b>	LEGIONELLA PNEUMOPHILA SG 1 ELISA IgG
<b>M1000</b>	LEGIONELLA PNEUMOPHILA SG 1 ELISA IgM
<b>T1001</b>	LEGIONELLA PNEUMOPHILA SG 1-6 ELISA IgG+IgM
<b>G1002</b>	MYCOPLASMA PNEUMONIAE ELISA IgG
<b>G/M1015</b>	MYCOPLASMA PNEUMONIAE ELISA IgG/IgM
<b>M1002</b>	MYCOPLASMA PNEUMONIAE ELISA IgM
<b>M1040</b>	PNEUMOBACT ELISA IgM <i>L. pneumophila</i> sg 1, <i>M. pneumoniae</i> , <i>C. burnetii</i> and <i>C. pneumoniae</i>

## VIRAL RESPIRATORY INFECTIONS

<b>G/M1004</b>	ADENOVIRUS ELISA IgG/IgM
<b>G/M1007</b>	INFLUENZA A ELISA IgG/IgM
<b>G/M1008</b>	INFLUENZA B ELISA IgG/IgM
<b>G/M1009</b>	PARAINFLUENZA 1 ELISA IgG/IgM
<b>G/M1010</b>	PARAINFLUENZA 2 ELISA IgG/IgM
<b>G/M1011</b>	PARAINFLUENZA 3 ELISA IgG/IgM
<b>G/M1006</b>	RESPIRATORY SYNCYTIAL VIRUS ELISA IgG/IgM

## TORCH INFECTIONS

<b>G1004</b>	CYTOMEGALOVIRUS ELISA IgG (includes calibrators)
<b>M1004</b>	CYTOMEGALOVIRUS ELISA IgM CAPTURE
<b>G/M1012</b>	HERPES SIMPLEX 1 ELISA IgG/IgM (gG1 recombinant)
<b>G/M1016</b>	HERPES SIMPLEX 1+2 ELISA IgG/IgM
<b>G/M1013</b>	HERPES SIMPLEX 2 ELISA IgG/IgM (gG2 purified)
<b>G1026</b>	RUBELLA ELISA IgG (includes calibrators)
<b>M1026</b>	RUBELLA ELISA IgM CAPTURE
<b>G1027</b>	TOXOPLASMA ELISA IgG (includes calibrators)
<b>M1027</b>	TOXOPLASMA ELISA IgM CAPTURE

## VACCINE-PREVENTABLE INFECTIONS

<b>G/M1001</b>	MEASLES ELISA IgG/IgM
<b>G/M1014</b>	MUMPS ELISA IgG/IgM
<b>G1008</b>	TETANUS ELISA IgG (includes calibrators)
<b>G/M1002</b>	VARICELLA-ZOSTER ELISA IgG/IgM

## VECTOR-BORNE INFECTIONS

<b>T1020</b>	CHAGAS ELISA IgG+IgM
<b>G1010</b>	CHIKUNGUNYA ELISA IgG
<b>M1010</b>	CHIKUNGUNYA ELISA IgM
<b>G1018</b>	DENGUE ELISA IgG
<b>M1018</b>	DENGUE ELISA IgM CAPTURE
<b>T1000</b>	LEISHMANIA ELISA IgG+IgM
<b>G/M1005</b>	RICKETTSIA CONORII ELISA IgG/IgM
<b>G1023</b>	ZIKA ELISA IgG
<b>M1023</b>	ZIKA ELISA IgM

## SEXUALLY TRANSMITTED INFECTIONS

<b>A1017</b>	CHLAMYDIA TRACHOMATIS ELISA IgA
<b>G/M1017</b>	CHLAMYDIA TRACHOMATIS ELISA IgG/IgM
<b>G1060</b>	SYPHILIS ELISA IgG
<b>T1060</b>	SYPHILIS ELISA IgG+IgM

## ZOONOTIC INFECTIONS

<b>G1003</b>	BRUCELLA ELISA IgG
<b>M1006</b>	BRUCELLA ELISA IgM
<b>G1006</b>	HYDATIDOSIS ELISA IgG

## OTHER INFECTIONS

<b>G1205</b>	EPSTEIN-BARR EA ELISA IgG
<b>G1105</b>	EPSTEIN-BARR EBNA ELISA IgG
<b>G1005</b>	EPSTEIN-BARR VCA ELISA IgG
<b>M1005</b>	EPSTEIN-BARR VCA ELISA IgM
<b>A1022</b>	HELICOBACTER PYLORI ELISA IgA (incl. calibrators)
<b>G1022</b>	HELICOBACTER PYLORI ELISA IgG (incl. calibrators)
<b>M1003</b>	LEPTOSPIRA ELISA IgM
<b>G1031</b>	PARVOVIRUS ELISA IgG
<b>M1031</b>	PARVOVIRUS ELISA IgM CAPTURE

## OTHER ELISA PRODUCTS

<b>S001</b>	ELISA SORBENT (2 x 1.2 ml)
<b>S001-100</b>	ELISA SORBENT (100 ml)
<b>B003</b>	SERUM DILUENT FOR ELISA KITS (500 ml)
<b>B001</b>	WASH BUFFER FOR ELISA KITS (500 ml 20x)

### Vircell, S.L.

Parque Tecnológico de la Salud, Avicena 8. 18016 Granada, Spain  
+34 958 441 264 | [info@vircell.com](mailto:info@vircell.com) | [www.vircell.com](http://www.vircell.com)